

# i-Vinci SCANNER

Rizhao | Beijing | Shanghai | Qingdao | LA

Manufactured for:

**VINCISMILE GROUP LLC**

Add: 20524 Carrey Rd. Walnut, CA 91789, USA  
Email: [marketing@vincismile.com](mailto:marketing@vincismile.com)  
Web: [www.vincismile.com](http://www.vincismile.com)

Manufacturer:

Ziyang Freqty Medical Equipment Co., Ltd.  
Add: Zone A, B, C, 4 Floor, Building A, No.3, Xiandai Road,  
Ziyang City, Sichuan Province, China

# i-Vinci Smart

Open A New Era of Digital Dentistry

Only **138g** 221x27x25mm



# SPECIFICATION

Technology	Staring Scanning with Coded Structure Light
Anti Fog System	Intelligent Heating
Accuracy	7 $\mu$ m
Precision	5 $\mu$ m
Handpiece Weight	138g (192g with cable)
Handpiece Dimension	221x27x25mm
Handpiece Enclosure	Aviation Aluminum Alloy
Scan Field	18x16mm
Scan Depth	0~20mm
True Color	Full HD
Tips	x6 ( Autoclave )
Tip Autoclave Cycle	50~100 times
Scanning Control	Button
Motion Sensing	Support
Data Interface	USB3.0
Cable Length (m)	1.8m

Smaller Tip

Faster Scan

Higher Accuracy

Motion Sensing

Super Quiet

